

electronic pipette controller

Profiller™ *electro* 446

Pipette controller fitting any hand in a very comfortable way. Shape is ideal for working with glass or plastic pipettes within laminar flow or microbiological safety cabinet. A powerful, yet quiet motor and a switch to set both aspiration and dispensing speeds enable maximal working efficiency even with large pipetting volumes. One-year warranty.

Product advantages

- Powerful pump - fills 25 mL in 3 seconds
- Simple, fast operation
- Comfortable finger triggers
- Low battery light
- NiMH battery rechargeable during use
- Interchangeable protection membrane filter
- Stand suitable for bench or wall mounting



*Pipettes
up to 100 mL*



Finger triggers



Variable speed



Battery charge

Convenient stand

Suitable for bench or wall mounting, stand houses instrument with or without a pipette in place.



Ordering information - Instrument

Description	Plug style	Pack.	Cat. No.
Profiller™ 446	Australia - N.Z.	1 / pk	446.100A
Profiller™ 446	Europe	1 / pk	446.100E
Profiller™ 446	United Kingdom	1 / pk	446.100G
Profiller™ 446	USA - Japan	1 / pk	446.100U

Supplied with charger, two spare hydrophobic filters (0.45µm and 0.2µm), bench/wall stand and instruction manual.

Ordering information - Accessories

Description	Pack.	Cat. No.
Spare bench stand for Profiller™ 446	1 / pk	320.446
Replacement 1.2V NiMH battery	1 / pk	900.916
Fast charger + cord, Europe style plug	1 / pk	900.946E
Fast charger + cord, UK style plug	1 / pk	900.946G
Fast charger + cord, USA-Japan style	1 / pk	900.946U
PTFE Membrane filter, 0.45 µm	5 / pk	322.435
PTFE Membrane filter, 0.2 µm	5 / pk	322.400
Spare silicone pipette holder	1 / pk	1.435.01
Spare nose cone white	3 / pk	1.435.02
Alt. nose cones pink, violet, green (1 ea.)	3 / pk	1.435.03

For more information or for an overview of the complete Socorex product range, visit www.socorex.com